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Biosint[®] Extra

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso

DeguDent
A Dentsply Company

Biosint[®] Extra

| | | |
|------|-------------------------------------|----|
| (D) | Gebrauchsanweisung | 5 |
| (GB) | Instructions for use | 11 |
| (F) | Mode d'emploi | 17 |
| (I) | Istruzioni per l'uso | 23 |
| (E) | Instrucciones para el uso | 29 |

Biosint® Extra

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Biosint® Extra

Biosint® Extra is a gypsum-free, phosphate-bonded special investment compound for the production of duplicate casts and dies for the CoCr model cast technique.

Biosint Extra is mixed with Biosol® E.

Adding the special Biosol E mixing liquid to water changes the total expansion.

Mixing only with water (distilled water) is possible.

1 Safety information

Mixing with other liquids or materials during the use must not be carried out, since during this process health hazards as well as impairments to the casting quality cannot be excluded.

- This product contains quartz and cristobalite, some in the form of fine dust that may penetrate the pulmonary alveoli. Inhalation of the dust constitutes a potential carcinogenic hazard. Ensure that an adequate suction line/ventilation is available in the workspace and near the machines.
- Inhalation the dust of the product may trigger silicosis.
- Employers must provide Class P3 respiratory protection.
- Heating the investment compound above 200 °C will result in the formation of ammonia. Ammonia irritates the respiratory organs.
- May cause mechanical irritation of the skin and mucosa of the eye and respiratory passages.
- Irritation and caustic reactions to the skin and mucosa must be expected on local contact.

Carry out the "coarse" removal of the casting after the casting operation under water or soak the investment material mold.

1.1 Information for the processing of mixing liquids for dental investment materials

Following contact with skin:

Upon skin contact, wash with plenty of water.

Following eye contact:

Rinse thoroughly with plenty of water keeping eyelid open.

In case of persistent discomfort: Consult an ophthalmologist.

Mixing liquids are alkaline.

The Safety Data Sheet includes further details.

2 Technical data

| | |
|--------------------------|--|
| Mixing ratio | 100 g powder : 15 ml liquid |
| Total expansion | 0.95 % – 1.65 % |
| Mixing time under vacuum | 60 sec |
| Processing time span | 3 – 5 min |
| Compression strength | according to concentration of the expansion liquid between 15 and 20 MPa |

The total expansion, i.e. the sum of setting expansion and thermal expansion can be changed by diluting the Biosol E mixing liquid.

| Mixing liquid consisting of: | | Concentration | Total sum of: | |
|------------------------------|----------|---------------|---------------|-------------|
| Water | Biosol E | Biosol E | setting exp. | therm. exp. |
| 4 parts | 0 parts | 0 % | 0.30 % | 0.65 % |
| 2 parts | 2 parts | 50 % | 0.35 % | 0.95 % |
| 1 parts | 3 parts | 75 % | 0.45 % | 1.05 % |
| 0 parts | 4 parts | 100 % | 0.50 % | 1.15 % |

3 Important information

The safety information mentioned under 1 and 1.1 must be observed prior to the usage.

Biosint Extra must not be brought into contact with gypsum or gypsum-containing materials since this may lead to sulphur damage of high-melting alloys. Therefore special mixing cups and spatulas must be used for mixing Biosint Extra.

3.1 Investment storage

Most of you are probably aware of the problems which can occur when investment is subjected to variations in environment temperature (e.g. summer and winter) and the effect this can have on your castings. That is why it is very important to find the right place to store your supply of investment. The investment should be stored in a cool and dry place, between +5 °C to +30 °C at 30 % to 90 % relative humidity.

If you have a place in your laboratory where these conditions can be found throughout the year, then that's the place to store your investment!

If it is not possible for you to store your investment at this temperature, then you should at least make sure that the work bench where investing is done fits these requirements (use a temperature control cupboard or similar).

Please note the symbols on product labels:

REF Product code

LOT Batch or lot number

 Use before

 Consult instructions for use

Avoid any nasty surprises by observing the following:

- Regularly check the temperature and humidity in the embedding region -**application temperature: +20 °C to +23 °C at 30 % to 90 % relative humidity.**
- Avoid using any investment that has just been brought out of an especially warm or cold place. (e.g. storage rooms, silos, cupboards etc.)!
- Fill up stocks at your workbench before your supply runs out. Give the investment time to adjust to the environment!

- Avoid using investment that has just been delivered to the laboratory. Again, give the investment time to acquire the right temperature.
- If space allows, stock up on investment (especially LIQUID) before the cold winter months arrive. This way, you forego the risk of exposing the liquid to frost during delivery

4 Production of duplicating models

4.1 Mixing of the investment compound

Premix the powder and liquid by hand with a spatula in the clean, freshly rinsed cup of the Multivac® evacuating unit. Then let the investment compound stand under vacuum for a short time. Subsequently mix under vacuum for approx. 60 sec.

4.2 Pouring the duplicating mould

Mix the investment compound flow into the duplicating mould under slight vibration and carefully avoid the formation of bubbles. After the duplicating mould is filled, immediately switch off the vibrator.

Note: Silicone duplicating mould can be poured easier if they have been wetted previously with Waxit®. The mould must be blown dry before pouring with investment compound.

Investment compound models which have to be produced in silicone duplicating moulds must be mixed with higher concentration of Biosol E; they require greater expansion.

4.3 Setting

Let the investment compound set for approx. 30 min and then remove carefully from the duplicating mould.

4.4 Drying

Dry the model for approx. 40 min in a drying furnace preheated to around 220 °C.

4.5 Hardening

Note: Investment compound models which have been produced in silicone duplicating moulds need not to be hardened.

Immerse the dried investment compound model in cold immersion hardener for approx. 10 s.

Let the immersion liquid drip off the cast and put the cast back in the switched-off drying furnace for approx. 5 min.

Using hardening waxes, hardening sprays or wax adhesives is also possible.

4.6 Modelling

Modelling on hand-warm model.

5 Production of the casting mould

5.1 Relaxing

Thinly brush the wax object with the Waxit wetting agent and then carefully blow it dry.

5.2 Fine investment

Before investing of the wax modelation with Biosint Extra, the entire wax surface can be pretreated with fine investment compound. This provides the best guarantee of avoiding air bubbles at undercuts or points on the poured object which are difficult to see and also results in exact surface reproduction of the cast.

5.3 Mixing of the investment compound

Dose Biosint Extra and mixing liquid and mix as described under 4.1.

We recommend as mixing ratio 100 g powder : 15 ml liquid.

Note: The casting mould should be produced in the same or somewhat lower concentration as the model!

5.4 Investing

Let the mixed investment compound flow into the mould form under light vibration and carefully avoid the formation of bubbles.

Immediately switch off the vibrator after the mould is filled.

5.5 Setting

Let the mould set for approx. 30 min.

5.6 Burn out of the wax

Place the mould in a furnace preheated to 300 °C and burn out the wax for approx. 45–60 min.

5.7 Preheating

Then bring the mould to the required final temperature of 1000–1050 °C required for CoCrMo alloys. Maintain the final temperature for approx. 60 min.

The notes provided under 5.6 and 5.7 also apply accordingly for program-controlled preheating overnight.

5.8 Casting

Cast according to the processing instructions of the alloy manufacturer.

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DeguDent GmbH